



UNO 5 Oxygen Concentrator

UNO Supply introduces the latest product in the Oxygen Concentrators category. UNO 5, an innovative and evolved product, is the result of our pursuit to provide the absolute best to our customers.

UNO 5 is lighter and more compact, with more functionality and benefits. Weighing only 15.2 kg, it is among the lightest in the industry and guarantees complete user satisfaction and convenience.

Technical Parameters

Delivery Rate (LPM)	0.5 - 5
Maximum Flow (LPM)	5
Output Pressure (KPa)	58 ± 6 (8.5 ± 0.5 psi)
Electrical Requirement Voltage Frequency	110VAC , 220VAC 50/60Hz
Purity (%)	93 ± 3 (0.5-5.0 L)
Power Consumption (Watts)	300
Weight (Kg)	15.2
Noise Level dB(A)	front 1 M < 42
Dimensions (LxWxH cm)	3.68 x 3.42 x 5.72
OCI	Optional

